

PFA Insulated IEC Fine Gauge Exposed Junction Thermocouple (0.076mm conductors) - Type K - 300mm Length





Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.





Thermocouple type K. Fast response, welded exposed junction available with bare tails. Ideal for test & development applications.

Specifications

Specifications

XE-0453-001
Fast response, welded tip exposed junction
Exposed Junction
K IEC
300mm
Teflon® PFA insulated wires
0.003" (0.076mm) conductors
Bare Tails
260°C
-75°C
IEC-584-3
Exposed Junctions