

PTFE Single Shot IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 300mm Cable Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Specifications

Specifications

e single shot cable construction Resistant to chemicals, oils and fluids. le for aerospace applications and is wildly used in autoclaves. led Junction
ed Junction
า
6mm
ure flat pin plug, colour coded in accordance IEC 584
14-3
ure
ed Junctions