

Glassfibre insulated Flat Pair IEC Exposed Junction Thermocouple - Type K with Miniature Plug Termination - 1/0.315mm x 300mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order, you will be contacted regarding delivery information

Specifications**Specifications**

Product Code	XE-0430-001
General Description	Good temperature resistance (not suitable in fluids)
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	300mm
Cable Type	Glassfibre
Core / Strands	1/0.315mm
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Colour	Green



Series	Exposed Junctions
--------	-------------------