

WT9055 Digital Sound Level Meter



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

The WT9055 digital sound level meter is novel, small and portable. The device is applicable for measurement of factory environments, quality control, health prevention and various environmental noise, including offices, transporting routes, stereo equipment and many more.

Sound level measurement(dBA)

MIN / MAX / Lock current value

Data Hold

LCD backlight function

Manual / auto shutoff

Backlight alarming

Software for these products can be downloaded by following

[this link](#)

Specifications

Specifications

Product Code	XE-4157-001
Model	WT9055
Type	Sound Level Meter
Accuracy	+/- 1.5dB
Resolution	0.1dB(A)
Battery	3 x 1.5V AAA
Dimensions	53.4mm x 31.4mm x 131.7mm
Working Temperature	0-40°C
Operating Humidity	10-80% RH
Weight	96.38g
Storage Temperature	10-60°C
Measuring range	35-130dB (A)



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Important Information

Frequency response: 31.5Hz - 8KHz