

## Micro Pocket Thermometer - White



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature instant-read thermometer

Compact design, perfect for chefs' pockets

Effortless transport during line checks

Reduced tip for quicker response

Temperature range: -50 to 150 °C

IP65 water-resistant

Designed for Speed and Convenience

This miniature instant-read thermometer is built for quick, convenient measurements on the go. Whether you're ensuring foods are hot for service or conducting food storage checks, the Micro Pocket Thermometer will always be by your side. At just 124 mm long, it will easily fit inside a pocket, making it the ideal companion for busy chefs who can't afford to waste time searching for their thermometer during peak hours.

HACCP Temperatures Covered

With a temperature range from -50 to 150 °C, this thermometer addresses all food safety needs, including chilled food storage and hot food service. Accurate to  $\pm 0.5$  °C (between -10 and 100 °C), it gives you confidence in your readings.

### Specifications

## Specifications

<b>Product Code</b>	XE-4202-001
<b>Accuracy</b>	$\pm 0.5$ °C (-10 to 100 °C), otherwise 1.0 °C
<b>Resolution</b>	0.1 °C
<b>Air Temp: Range</b>	-50 °C to 150 °C
<b>Battery</b>	1 x 3 volt CR1632
<b>Colour</b>	White
<b>Dimensions</b>	22 x 27 x 124 mm
<b>Display Type</b>	13 x 9mm LCD
<b>Automatic Shutdown</b>	Auto off after 5 minutes
<b>Power Consumption</b>	2000 hours battery life
<b>Standards Met</b>	IP65



<b>Construction</b>	ABS plastic includes Biomaster product protection
---------------------	---