

ThermaPod USB temperature data logger



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Rapid software-free data downloads

Ideal for storage and transportation

Records up to 7,500 readings

Splashproof to IP65

Temperature range -20 to 85 °C
 Automated Temperature Recording

Automate ambient temperature monitoring affordably with the ThermoPod USB logger. With a range of -20 to 85 °C, it's ideal for monitoring storage conditions for chilled items like food and pharmaceuticals or ambient temperatures in rooms, storerooms, and zoos.

With no software download is required, it's perfect for processes like transportation. Users can generate a PDF report instantly on arrival, allowing quick confirmation on whether goods should be accepted or rejected.

Easy to Set Up

Program your set intervals easily with two buttons on either side of the display. Choose from 6 to 255 minutes. The LED flashes green when recording. The logger can store up to 7,500 readings.

Specifications

Specifications

Product Code	XE-4204-001
Accuracy	+/-0.5°C (-10 to 70°C) otherwise +/-1°C
Readings	7,500 memory
Resolution	0.1°C
Logging Interval	6 seconds to 255 minutes
Battery	3 volt CR2032 lithium coin cell
Colour	White
Dimensions	14.5 x 26 x 76 mm
Automatic Shutdown	Auto off after 30 seconds
Max. Temperature	85°C
Min. Temperature	-25°C
Weight	26 grams



Power Consumption	1000 hours battery life
Standards Met	IP65