

## RFX Wireless Meat Probe



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Track core meat temperatures from your device

Use up to 50 probes with one GATEWAY

Reliable wireless connection

50 hours of battery life with a 10-minute charge

Fully waterproof to IP69K

Free traceable calibration certificate

*Must be paired with the*  
*or*

[RFX Gateway](#)  
[RFX Starter Kit](#)

### Specifications

## Specifications

Product Code	XE-4201-001
Model	RFX Wireless Meat Probe
Type	Thermistor, 4 Sensors
Accuracy	+/-0.5°C (-10 to 100°C)
Battery	1.5 V AAA with charging block - 175+ full charges per AAA battery. 0-80% in 10 minutes; 0-100% in 20 minutes
warranty	6 Months
Output	Wireless communication: 433Mhz radio frequency
Temperature Range	Internal: -18 to 100°C / External: -18 to 538°C
Standards Met	IP69K
Measuring range	Battery powered RFX GATEWAY: 650m (2,132 ft.), AC Powered RFX GATEWAY: 385m (1,263 ft.) (RFX MEAT to RFX GATEWAY)