



Miniature Rectangular Fascia Socket Thermocouple Connector AM-K-FF Type K ANSI



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Rectangular Thermocouple Connector Fascia Socket AM-K-FF Type K ANSI

Nylon clip for quick, simple fixing to panel. All fascia mounting sockets are polarised by means of unequally sized pins

Specifications**Specifications**

Product Code	XE-1137-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	K ANSI
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Yellow
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector