



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI K Needle Tip Temperature Probe - 3.3 x 110mm with Mini Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

This needle tip Temperature probe (3.3 × 110 mm) with mini plug is designed for fast and accurate temperature measurement in soft and semi-solid materials. The fine stainless-steel needle tip provides a rapid response and allows easy insertion with minimal disturbance to the sample. Durable and versatile, it's suitable for food testing, laboratory work, and a wide range of general-purpose applications.

Specifications

Specifications

Product Code	XE-0043-001
General Description	Needle Tip
Thermocouple Type	ANSI
Probe Diameter	3.3mm
Probe Length	110mm
Termination Type	Miniature Plug
Max. Temperature	260°C
Min. Temperature	-50°C
Standards Met	ANSI Type K
Connector Type	Plug
Connector Size	Miniature
Colour	Yellow