




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI K Air Tip Temperature Probe 4 x 110mm with Mini Plug



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

This air tip temperature probe (4 × 110 mm) with mini plug is designed for rapid and accurate measurement of air and gas temperatures. Its perforated or exposed-tip design allows excellent airflow around the sensor, delivering fast response times and stable readings. Robust and versatile, it's ideal for HVAC testing, environmental monitoring, and a wide range of laboratory or industrial applications.

Specifications

Specifications

Product Code	XE-0042-001
General Description	Air Tip Probe
Thermocouple Type	ANSI
Probe Diameter	4mm
Probe Length	110mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Standards Met	ANSI Type K
Connector Type	Plug
Connector Size	Miniature
Colour	Yellow