


SCH4/ABS Terminal Head + 2 x M16 glands



M16 Nylon Cable Gland

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

SCH4/ABS Terminal Heads provide durable, impact-resistant, and corrosion-resistant enclosures for protecting electrical connections in various applications. The robust ABS construction, combined with the practical design and compatibility, makes them suitable for a wide range of industrial, commercial, and HVAC applications.

Proper installation and regular maintenance of SCH4/ABS Terminal Heads ensure reliable performance and long-term protection of sensitive electrical components. The SCH4/ABS head has 2 x M16 entry threads.

Specifications

Specifications

Product Code	XE-3842-001
Head Type	SCH4/ABS Terminal Head
Material	ABS Plastic
Cable Entry	M16 (x2)
Fixing (Centres)	33mm & 40mm
Max. Temperature	70C
IP Rating	IP67 rated
Cap Retaining Type	Screw cap and retaining strap