


Miniature Panel Mount Thermocouple Connector with Stainless Steel Bracket AM-K-SSPF Type K ANSI



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Panel Mount Thermocouple Connector with Stainless Steel Bracket AM-K-SSPF Type K ANSI

Stainless steel bracket simply fitted to connector. Fully enclosed terminals with cable port seal

Specifications**Specifications**

Product Code	XE-1138-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	K ANSI
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector