

Lascar EL-PROBE-EXTENDERS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

EL-PROBE-EXTENDERS

3.5mm jack to female 3.5mm socket which can be used to extend the probe supplied with the all of Lascar's Thermistor (TP and TP+) products.

Specifications**Specifications**

General Description	Probe extenders for Lascar data loggers
Type	3.5mm jack to female 3.5mm socket
Cable Length	1.5M, 3M, 5M, 10M