

BF Wireless Meat Thermometer with Wi-Fi & Bluetooth



Cook With Confidence, Not Constraint



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Single Probe (BF70) or Dual Probe (BF80) options. Waterproof and Heat-Resistant.

This BF80 wireless meat thermometer has two separate stainless steel probes to measure two foods simultaneously and improve BBQ cooking efficiency. Its heat-resistant probes adapt to ovens and grills, and you can rinse them with cool water after they cool completely post-cooking. BF70 version has the same features with a single probe.

Dual Wireless Connection. As a multi-functional BBQ thermometer probe, it is built with dual wireless modules: 2.4GHz Wi-Fi enables boundless remote oversight no matter how far away you are from cookers, whereas upgraded Bluetooth 5.4 offers steady signal coverage up to 150 meters. This handy meat probe thermometer keeps you informed of cooking status flexibly indoors or outdoors.

150° Wide View Screen & Precise TI Chip. Fitted with a 150° panoramic wide-view positive display and high-performance TI chip, this wireless meat thermometer smart base delivers ultra-stable, accurate temperature readings for grilling and baking. It splits the screen to display ambient heat and pre-set temperature, supporting simple offline touch operation without a phone.

Custom APP Management. You can fully customise various meat target temperatures through the exclusive matching APP of this BF wireless meat thermometer. The system automatically generates visualised temperature tracking maps to save your cooking records, helping you tweak cooking skills and optimise homemade barbecue recipes step by step.

Long Battery & Versatile Gift. The base offers 24-hour runtime, and probes get 33 hours of use after a 10-minute quick charge, enough for all-day low-temp smoking and lengthy home roasting. Compatible with ovens, air fryers and outdoor grilling scenarios, this wireless meat thermometer turns into a practical surprise gift for grill fanatics and household cooking lovers on holidays.

Specifications

Specifications

Model	BF70 (single probe) or BF80 (Dual probe)
General Description	Wireless Meat Thermometer
Communications Protocol	WiFi & 100m Bluetooth Range
Accuracy	+/- 1°C (+/- 1.8F)
Ambient Operating Temperature	up to 400°C (752°F)
Battery	1 x 9V battery (included). Rechargeable
Display Type	2.4 inch LCD display
Power Consumption	Wireless probes: 5 min fast charging for 33+ hours of use