



Standard Rectangular Fascia Thermocouple Connector Socket JS-R/S-FF Type R/S JIS



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Rectangular Thermocouple Connector Fascia Socket JS-R/S-FF Type R/S JIS

Nylon clip for quick, simple fixing to panel. All fascia mounting sockets are polarised by means of unequally sized pins

Specifications**Specifications**

Product Code	XE-9317-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	R/S JIS
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Black
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS