


## Strain Relief Grommet



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Strain Relief Grommet for Miniature Connectors**

Temperature Rating: -30°C to +140°C

Material: Thermoplastic elastomer (TPE)

**Specifications**

<b>Product Code</b>	XE-1882-001
<b>Colour</b>	Black
<b>Max. Temperature</b>	+220°C
<b>Min. Temperature</b>	-35°C
<b>Series</b>	GROMMET